


R2GATE™

turning imagination into reality

since 2012



7 R2GATE™ guided surgery

7.1 Guided surgery

Pre-surgery checkpoints

1
Pre-surgery
checkpoints

R2 Guide
fitting

Drilling
protocol

ONE-DAY
IMPLANT

R2
GATE

Full Manual
ver20180328

R2 Guide™

Virtual planning to reality!

- 3D printed based on approved treatment plan
- One body combining drill guide hole(s), drill stopper(s), & hex control



R2 guided surgery

Highly accurate & convenient

- All drills include drill, drill guide & drill stopper as one-body
- No need for metal sleeves or spoons!
- Shorter surgery time!



1 Sterilization of R2 Guide & prosthesis

Sterilize R2 Guide & prosthesis in **ANTISEPTIC** (e.g. chlorhexidine or gluconate) for 30 minutes before surgery



! Do NOT use autoclave

- Autoclave causes distortion of R2 Guide
- Sterilize R2 Guide using antiseptic



2 Wide implant (over Ø5.0mm - real size)

In this case, R2 Centre provides two R2 Guides:

- regular drill guide hole (Ø5.0)
- wide drill guide hole (Ø5.0-6.5)

Drill to 4.8mm drill using regular drill guide hole, then use wide drill guide hole to finish drilling & place wide implant



7.2 Guided surgery

R2 Guide fitting



Pre-surgery checkpoints



R2 Guide fitting



Drilling protocol





ONE-DAY IMPLANT

1 Tooth & tissue support

For improved accuracy, R2 Guide is seated using several seating holes for neighboring teeth

Once R2 Guide is seated, check fit of seating holes & cusp tip of neighboring teeth



2 Tissue-only support (edentulous case)

In this case, R2Guide is seated using putty bite & specially designed anchor pins



Anchor Kit:

R2 Guide is fixed using anchor pins inserted via pinholes in R2 Guide

Triangular format of anchor pins is recommended for better fixation.



Place R2 Guide & putty bite in mouth together



Bite with maximum occlusal force



Screw anchor pins into pin holes using hand driver



Advance drilling (2.0 x 13mm) can be used in case of hard bone (mandible)

7.3 Guided surgery Drilling protocol



Pre-surgery
checkpoints



R2 Guide
fitting



Drilling
protocol



ONE-DAY
IMPLANT

**R2
GATE**

Full Manual
ver20180328

R2GATE™

Compatible with all major implant systems



R2 Surgical Kits

Various kits available!



1 R2 Full Kit

Complete set of drills + system-specific implant carriers for MegaGen implant systems [AnyRidge, AnyONE]



2 R2 Universal Kit


Essential drills (initial - Ø2.8 drills) for all major implant systems
Implant carriers & final drills can be added for preferred implant system



3 R2 Narrow Kit

For MINI implant system:
Ø3.5mm drill guide is designed for narrow surgical sites, such as anterior mandible or close adjacent teeth

7.3 Guided surgery Drilling protocol




R2 Guide drilling protocol

- All drills combine drill, drill guide, and drill stopper in one-body
- No need for metal sleeves or spoons!
- Shorter surgery time!
- Disposable final drill provided for each surgery to optimize initial stability



1 Narrow crest drill for narrow or steep alveolar ridge

To prevent slipping in cases of narrow or steep alveolar ridge, use narrow crest drill first to flatten drilling area



- a** Engage bone when drilling counter-clockwise at less than 100 RPM
- b** Drill clockwise at 400~600 RPM until drill is stopped by drill stopper
- c** Withdraw drill
- d** Mark drilling start point on bone using initial twist drill
- !** Autogenous bone can be harvested from shank of narrow crest drill




R2 Universal Kit



R2 Guide drilling protocol


2 Initial & 2nd drilling

- Use initial drill to mark drilling start point on bone
- 2nd drill also removes excess bone above fixture platform for better connection with prosthesis
- In case of dense bone, use 2nd drill right before fixture placement




3 Graduated Drilling: **critical to success!**

- Use graduated drilling sequence to increase osteotomy depth & diameter
- E.g. when placing 11.5mm fixture: narrow drill -> initial drill -> 2nd drill -> 2.0x7 -> 2.0x8.5 -> 2.0x10 -> 2.0x11.5 -> 2.5x11.5 -> 2.8x11.5 -> final drill -> cortical bone drill




4 Drill guide: **for maximum accuracy**

- Do NOT start drilling until drill guide is completely inserted into guide hole
- Always start drilling at low RPM [300~500 RPM]




5 Slow UP & DOWN drill motion

Increase osteotomy depth in 1mm increments using slow up & down motions until drill stopper reaches stopper position in drill guide



7.3 Guided surgery Drilling protocol



R2 Guide drilling protocol



6 Drilling strategy according to bone density

Color-coded analysis of bone morphology identifies invisible bony structure for predicting optimal drilling sequence to achieve strong initial stability

Recommended drilling speed **500 - 800 RPM**

Symbols used in drilling protocol report

- Full-depth drilling
- ◐ Half-depth drilling
- !! Optional drilling depending on decision of clinician



7.3 Guided surgery Fixture delivery



Pre-surgery
checkpoints



R2 Guide
fitting



Drilling
protocol



ONE-DAY
IMPLANT


R2
GATE

Full Manual
ver20180328


7 Fixture delivery using implant carrier

Fixture is delivered via R2 Guide using
implant carrier

- a** Fixture depth control:
Align upper line of implant carrier with R2
Guide window




- b** Hex position control:
Align green part of implant carrier with R2
Guide window turning in buccal direction



Implant placement in fully edentulous case

Anchor screws

In addition to anchor pins, use of
anchor screws in triangular
formation is highly recommended
in fully edentulous case to
increase stability of R2 Guide



Regular & wide fixtures

- When placing regular & wide
fixtures, R2 Centre will provide
2 R2 Guides
- Place regular fixtures first & use
same anchor screw positions
for placing wide fixtures



Anchor screws

7.4 R2 Guide Surgery ONE-DAY IMPLANT



ONE-DAY IMPLANT™

An effective reality!

R2GATE treatment planning enables pre-fabrication of customized drilling guide & prosthetics before surgery, which reduces chair-time, minimizes surgical procedures & allows for immediate loading

Customized prosthetics

Customized abutment & crown



- ZrGEN (Ti-Base)
- customized abutment
- zirconia/titanium provisional crown

Screw-retained crown (All-in-One)



- ZrGEN (Ti-Base)
- PMMA CAD/CAM provisional crown



All-on-Four(Six) type



- customized abutment
- multi-unit abutment
- 3D printed or milled dentures



Recommended conditions for immediate loading (ONE-DAY IMPLANT)

Based on extensive clinical trials, two key values need to be checked [limited to AnyRidge system]:

a Insertion Torque Value (ITV): >45Ncm

b Implant Stability Quotient (ISQ): > 70



! If either value is insufficient, loading should be delayed

ISQ & ITV conditions for various restorations

