

Immediate Denture – Procedure

In Clinic


In Lab

• Send Impression, Opposing, Bite to Lab

• **Make Record Base Plate (Resin Base),
Wax Bite Block (Wax Rim) : 3 days**

• Bite Registration
- Use Record Base Plate, Wax Bite Block
• Midline / Lip Line - Decision
• Shade Decision

• We can finish here, if doctor want.
• If doctor wants to do Try-In, do below procedure:
➤ Bite Check
➤ Teeth Setup : 3 days
➤ - Teeth Type
➤ - Shade match




• Finish : 3 - 5 days

Not Fit

Fit


Full Denture – Procedure



- *Shade Guide - Type*
 - Vita : A, B, C, D series (Recommended)
 - Bioform : B59, B62, B66, (Recommended)
 - Kenson : 55, 58, 63, (NOT recommended)

- *Teeth Type*
 - Regular (Basic) : Dentsply Classic, Yamahachi Naperce
 - Premium : Dentsply PN, iVoclar
 - Porcelain Teeth : Vita, iVoclar


Partial Denture – Procedure



- *Shade Guide - Type*
 - Vita : A, B, C, D series (Recommended)
 - Bioform : B59, B62, B66, (Recommended)
 - Kenson : 55, 58, 63, (NOT recommended)

- *Teeth Type*
 - Regular (Basic) : Dentsply Classic, Yamahachi Naperce
 - Premium : Dentsply PN, iVoclar
 - Porcelain Teeth : Vita, iVoclar

Locator OverDenture – Procedure



Needed Parts (*)
 - Fixture Level Impression: Impression Coping, Analog, (Locator Abmt, Male Processing Kit → Lab will decide)
~~- Abutment Level Impression: Locator Impression Coping, Locator Analog, Male Processing Kit~~

➤ Teeth Type
 - Regular (Basic) : Dentsply Classic, Yamahachi Naperce
 - Premium : Dentsply PN, iVoclar
 - Porcelain Teeth : Vita, iVoclar

➤ Shade Guide - Type
 - Vita : A, B, C, D series (Recommended)
 - Bioform : B59, B62, B66, (Recommended)
 - Kenson : 55, 58, 63, (NOT recommended)